

Experiment No. 3	Date:16 63 2021
[Module 3] Determination of FEC102.3: Interpret the basic known understanding the working of semiconduction [Module 2] FEC102.2: Explain the basic concepts of	quantum mechanics,
AIM: - To draw the V-I characteristic for I determine the value of Planck's constant.	
Apparatus: - Planck's constant kit, conne	cting wires.
Theory: - Planck's constant is the fundamelates the energy of a photon to its frequency light Emitting Diodes (LED). Diodes today is achieved by having a slightly different so the experiment is based on the fact that frequency as: E = hv Where, E = energy of photon, h = Plance of the emitted photon	ency. To determine this constant we use come in a variety of colours. Each colour semiconductor material. the energy of the photon relates to its anck's Constant,
The electric energy needed to switch on a Where, $e =$ electric charge of an electron $V_o =$ Threshold voltage required to junction. This energy is different for different LE v). Each electron-hole radiative recombinate $E = hv = h(c/\lambda)$ Where, $c =$ speed of light (3 x 10 ⁸ m/s) $\lambda =$ wavelength of light.	(e = -1.6 x10 ⁻¹⁹ Coulomb) covercome the barrier created by p-n D colour (frequency of emitted light,
This threshold voltage (V_0) is related to the E_g = e V_0 + ΔE	om one LED to another.
Eq. (1) becomes $hv = e V_o + \Delta E$	
h = Planck's Constant e =	electronic charge velocity of light

In this experiment we use Eq. (3) to determine Planck's constant h. We use several LED's that emit light of different colours (frequencies).

Circuit Diagram:-



Procedure:-

- 1. Take the Planck's constant Determination Setup. Make the connections as shown in circuit diagram.
- 2. Connect + ve terminal of power supply to + ve terminal of DC voltmeter and -ve terminal to ve terminal of DC voltmeter.
- 3. Now take another patch cord and connect + ve of power supply to + ve of red colour LED.
- 4. Connect ve of power supply to ve of ammeter.
- 5. Connect +ve of ammeter to -ve of selected LED.
- 6. Set the voltmeter at the range of 20 V and ammeter at the 200 mA.
- 7. Connect the mains cord and switch ON the power supply.
- 8. Now increase the DC voltage at the fix interval as given in the observation table.
- 9. Note the corresponding current by DC ammeter in observation table.
- 10. Now take the current on Y-axis and voltage on X-axis and plot a graph between current and voltage
- 11. Find the knee- voltage or threshold-voltage from the graph.
- 12. Now switch OFF the DC power supply and break the LED connection.
- 13. Repeat above experiment for different colours of LEDs.
- 14. Put this value in given formula and calculate the Planck's constant $h = eV_0\lambda/c$

Take mean value of h calculated for different LEDs.

Observation

Table I: I-V Characteristics of different colour LEDs

	1	Red	Orange	Yellow	Green	Blue
	Voltage	Current	Current	Current (mA)	Current (mA)	Current (mA)
Sr.No.	(V)	(mA)	(mA)	-	0	0
01	0	0	0	0	0	0
02	0.5	0	0	6	0	6
03	1	0	0	0	0	0
04	1-2	0	0	0	0	0
05	1-3	0	0	0	0	0
06	1.4	0	0	0	0	0
07	1.5	0	0	0	0	0
08	1.6	0	0	0	0	0
09	1.7	0.1	0.1	0.3	0	0
10	1.8	0.9	0.4	1.4	0	0
11	1.9	5.1	2.4	3.9	0	0
12	2.0	13.6	5.5	7.1	0	0
13	2.1	25.4	16.6	10.8	0	0
14	2.2	40.1	(6.6	14.1	0	0
15	2.3	-	-	17.7	0	0
16	2.4	-	-	23.2	0	0
17	2-5	-		26.5	0.1	0.1
18	2.6	-	-	30.9	0.5	0.8
19	2.7	-	-	_	1.5	2.6
20	2.8	-	-	-	2.9	4.3
21	2.9	-	-	-	4.6	6.3
22	3.0					

Table II: Determination of Planck's constant Electronic Charge $e = 1.6 \times 10^{-19}$ coulomb Velocity of Light $c = 3 \times 10^{8}$ m/s

LED	Threshold Voltage Vo (in V) from	Wavelen gth (Å)	$h = eV_0\lambda/c$ (in J s)	Mean h (in J s)
	graph 2.65	4700	6.643 ×10-34	
Blue	2.68	5250	7.504 ×10-34	6 . 5 606 XIO 34
Green	1.85	5800		B . 3 600 Me
Yellow	1.86	6300	6.25 × 10-34	
Orange Red	1.79	7000	6.683 ×10-34	

Result:

Standard Value of Planck's constant = 6.626×16^{34} IS

Calculated Value of Planck's constant = 6.5606 × 10 IS

COMMENTS:

1. What is the difference between normal PN Junction diode and LED?

TED emits light while PN Junction diode connot emit light.

-> In LED all energy is converted into light cohile in PN Junction diode all energy is converted into heat.

2. Give examples of semiconductor materials which are used as LED. breakdown voltage.

i GaAs (Gallium arsenide)

2) Gap (challium Phosphide)

D.J.S.C.E. (Phy	ysics)
Journal	
Knowledge	3
Documentation	3
Punctuality	3
Virtual Lab (Performance & Documentation)	6
Total	15

Date	Signature of the faculty

